

Anti-Scd5 Polyclonal Antibody For Wb, If/Icc, Elisa - Q86Sk9

Item Number: CM0013893

Introduction

Rabbit polyclonal antibody against human SCD5 (Stearoyl-CoA desaturase 5, alias SCD4/ACOD4) validated for WB, IF/ICC, and ELISA. Involved in fatty acid desaturation, neuronal proliferation, and deafness.

[Learn More](#)

Parameter	Value
Product Name	SCD5 Rabbit pAb
Remarks/Alias	SCD4; ACOD4; FADS4; HSCD5; DFNA79; SCD5
Species	Human
Gene ID	79966
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-75 of human SCD5 (NP_001032671.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q86SK9
Protein Weight	38kDa
Shipping	Ice bag