

Anti-Scd5 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Sk9

Item Number: CM0026609

Introduction

High-quality rabbit anti-SCD5 polyclonal antibody validated for Western blot and ELISA. Detects human, mouse, and rat SCD5, a key enzyme in fatty acid biosynthesis and neuronal differentiation.

[Learn More](#)

Parameter	Value
Product Name	SCD5 Rabbit pAb
Remarks/Alias	SCD4; ACOD4; FADS4; HSCD5; DFNA79; SCD5
Species	Human
Gene ID (Human)	79966
Gene ID	79966
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 231-330 of human SCD5 (NP_001032671.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86SK9
Protein Weight	38kDa
Shipping	Ice bag