

Anti-Scfd1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wvm8

Item Number: CM0025192

Introduction

Rabbit polyclonal antibody against SCFD1 (also known as SLY1, STXBP1L2) for Western blot, IHC-P, and ELISA. Validated in human, mouse, and rat. Ideal for studying ER-Golgi transport.

[Learn More](#)

Parameter	Value
Product Name	SCFD1 Rabbit pAb
Remarks/Alias	SLY1; RA410; SLY1P; STXBP1L2; C14orf163; SCFD1
Species	Human
Gene ID (Human)	23256
Gene ID	23256
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 373-642 of human SCFD1 (NP_057190.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WVM8
Protein Weight	72kDa
Shipping	Ice bag