

Anti-Scg2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P13521

Item Number: CM0027194

Introduction

Rabbit polyclonal antibody against human SCG2 (Secretogranin-2), validated for WB, IF-P, ELISA. Reacts with human and rat. Suitable for neuroendocrine research and secretory granule biogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	SCG2 Rabbit pAb
Remarks/Alias	SN, CHGC, EM66, SgII, SCG2
Species	Human
Gene ID	7857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-320 of human SCG2 (NP_003460.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P13521
Protein Weight	71kDa
Shipping	Ice bag