

Anti-Scgb1A1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P11684

Item Number: CM0023561



Introduction

SCGB1A1 Rabbit pAb detects human SCGB1A1 (Uteroglobin/CC10), a club cell secretory protein, validated for WB, IF-P, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for respiratory research.

[Learn More](#)

Parameter	Value
Product Name	SCGB1A1 Rabbit pAb
Remarks/Alias	UGB, UP1, CC10, CC16, CCSP, UP-1, CCPBP, SCGB1A1
Species	Human
GeneID (Human)	7356
GeneID	7356
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-91 of human SCGB1A1 (NP_003348.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P11684
Protein Weight	10kDa
Shipping	Ice bag