

Anti-Scgb3A2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96PL1

Item Number: CM0016541

Introduction

Rabbit polyclonal antibody targeting human SCGB3A2, a secretoglobin with anti-inflammatory and anti-fibrotic roles. Validated for WB, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Ideal for lung development and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	SCGB3A2 Rabbit pAb
Remarks/Alias	LU103, PNSP1, UGRP1, pnSP-1, SCGB3A2
Species	Human
GeneID (Human)	117156
GeneID	117156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human SCGB3A2 (NP_473364.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PL1
Protein Weight	10kDa
Shipping	Ice bag