


Anti-Scin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6U3

Item Number: CM0019908



Introduction

Rabbit polyclonal antibody against human SCIN (Scinderin/Adseverin) for WB, IF/ICC, and ELISA. Reacts with human and mouse. Suitable for exocytosis, megakaryocyte, and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	SCIN Rabbit pAb
Remarks/Alias	SCIN; adseverin
Species	Human
Gene ID	85477
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SCIN (NP_001106177.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6U3
Protein Weight	80kDa
Shipping	Ice bag