

Anti-Scly Polyclonal Antibody For Wb, Elisa - Q96I15

Item Number: CM0014523



Introduction

High-quality Anti-SCLY rabbit polyclonal antibody for Western blot and ELISA applications. Reacts with human, mouse, and rat SCLY. Suitable for research on selenocysteine metabolism.

[Learn More](#)

Parameter	Value
Product Name	SCLY Rabbit pAb
Remarks/Alias	SCL; hSCL; SCLY
Species	Human
GeneID (Human)	51540
GeneID	51540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 284-453 of human SCLY (NP_057594.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96I15
Protein Weight	48kDa
Shipping	Ice bag