

## Anti-Scml1 Polyclonal Antibody For Wb, Elisa - Q9Un30

Item Number: CM0019173

### Introduction

Rabbit polyclonal antibody against human SCML1, a putative Polycomb group protein involved in transcriptional repression. Validated for WB and ELISA, cross-reacts with rat. Ideal for spermatogenesis and developmental biology research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | SCML1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SCML1                                                                                                              |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 6322                                                                                                               |
| GeneID           | 6322                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human SCML1 (NP_006737.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Rat                                                                                                                |
| SWISS            | Q9UN30                                                                                                             |
| Protein Weight   | 37kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |