


Anti-Scn1B Rabbit Polyclonal Antibody For Wb, Elisa - Q07699

Item Number: CM0014446



Introduction

Anti-SCN1B rabbit polyclonal antibody validated for Western blot and ELISA. Recognizes human, mouse, and rat SCN1B, a voltage-gated sodium channel regulatory subunit involved in neuronal excitability and cardiac function. Ideal for neuroscience and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	SCN1B Rabbit pAb
Remarks/Alias	DEE52; ATFB13; BRGDA5; EIEE52; GEFSP1; SCN1B
Species	Human
Gene ID (Human)	6324
Gene ID	6324
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-160 of human SCN1B (NP_950238.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07699
Protein Weight	25kDa
Shipping	Ice bag