


Anti-Scn3B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ny72

Item Number: CM0014012



Introduction

Rabbit polyclonal antibody targeting human SCN3B, a sodium channel regulatory beta subunit. Validated for WB, IF/ICC, ELISA, and cross-reacts with mouse and rat. Useful for cardiac electrophysiology and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	SCN3B Rabbit pAb
Remarks/Alias	SCNB3; ATFB16; BRGDA7; HSA243396; SCN3B
Species	Human
GeneID (Human)	55800
GeneID	55800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-160 of human SCN3B (NP_060870.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NY72
Protein Weight	25kDa
Shipping	Ice bag