


Anti-Scn5A Polyclonal Antibody For Wb, Elisa - Q14524

Item Number: CM0019271



Introduction

High-quality rabbit polyclonal antibody targeting SCN5A (Nav1.5) sodium channel, validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for cardiac electrophysiology and channelopathy research.

[Learn More](#)

Parameter	Value
Product Name	SCN5A Rabbit pAb
Remarks/Alias	HB1; HB2; HH1; IVF; VF1; HBBB; ICCD; LQT3; SSS1; CDCD2; CMD1E; CMPD2; PFHB1; Nav1.5; SCN5A
Species	Human
Gene ID (Human)	6331
Immunogen	Synthetic peptide: a synthetic peptide corresponding to a sequence within amino acids 1950-2016 of human SCN5A (NP_932173.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q14524
Protein Weight	227 kDa
Shipping	Ice bag