

Anti-Scnn1B Polyclonal Antibody For Wb And Elisa - P51168

Item Number: CM0022106

Introduction

Rabbit polyclonal antibody targeting the epithelial sodium channel subunit beta (SCNN1B/ENaC beta) for WB and ELISA applications. Specific to human and mouse SCNN1B, ideal for renal and respiratory research.

[Learn More](#)

Parameter	Value
Product Name	SCNN1B Rabbit pAb
Remarks/Alias	BESC1, ENaCb, SCNEB, LIDLS1, PHA1B2, ENaCbeta, beta-ENaC, beta-NaCH, SCNN1B
Species	Human
Gene ID (Human)	6338
Gene ID	6338
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 406-640 of human SCNN1B (NP_000327.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51168
Protein Weight	73kDa
Shipping	Ice bag