

Anti-Scnn1G Polyclonal Antibody For Wb, Ihc-P, Elisa - P51170

Item Number: CM0024716



Introduction

Rabbit polyclonal antibody against human SCNN1G (ENaC gamma subunit). Validated for Western blot, immunohistochemistry, and ELISA. Cross-reacts with mouse. Ideal for studying sodium channel physiology, kidney function, and ion transport.

[Learn More](#)

Parameter	Value
Product Name	SCNN1G Rabbit pAb
Remarks/Alias	PHA1; BESC3; ENaCg; LDL52; SCNEG; PHA1B3; ENaCgamma; SCNN1G
Species	Human
Gene ID (Human)	6340
Gene ID	6340
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51170
Protein Weight	74kDa
Shipping	Ice bag