

Anti-Scnn1G Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P51170

Item Number: CM0020170



Introduction

Rabbit polyclonal antibody against human SCNN1G (ENaC gamma subunit). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for sodium channel and epithelial transport research.

[Learn More](#)

Parameter	Value
Product Name	SCNN1G Rabbit pAb
Remarks/Alias	PHA1; BESC3; ENaCg; LDLS2; SCNEG; PHA1B3; ENaCgamma; SCNN1G
Species	Human
GeneID	6340 (human)
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-250 of human SCNN1G (NP_001030.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51170
Protein Weight	74kDa
Shipping	Ice bag