


Anti-Sco1 Polyclonal Antibody For Wb, Ihc-P, Elisa - 075880

Item Number: CM0025511



Introduction

Rabbit polyclonal antibody against human SCO1, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for mitochondrial copper metabolism and cytochrome c oxidase assembly studies.

[Learn More](#)

Parameter	Value
Product Name	SCO1 Rabbit pAb
Remarks/Alias	SCOD1; MC4DN4; SCO1
Species	Human
GeneID (Human)	6341
GeneID	6341
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human SCO1 (NP_004580.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075880
Protein Weight	34kDa
Shipping	Ice bag