

## Anti-Sco2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43819

Item Number: CM0029501



### Introduction

Rabbit polyclonal antibody against human SCO2, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat SCO2 (30 kDa), a mitochondrial copper chaperone essential for cytochrome c oxidase assembly and respiratory chain function.

[Learn More](#)

| Parameter        | Specification                                                                     |
|------------------|-----------------------------------------------------------------------------------|
| Product Name     | SCO2 Rabbit pAb                                                                   |
| Remarks/Alias    | TP; MYP6; TYMP; ECGF1; SCO1L; MC4DN2; CEMCOX1; PD-ECGF; TdRPase; Gliostatin; SCO2 |
| Species          | Human                                                                             |
| GeneID (Human)   | 9997                                                                              |
| GeneID           | 9997                                                                              |
| Immunogen        | Recombinant protein                                                               |
| Source           | Rabbit                                                                            |
| Category         | Polyclonal Antibodies                                                             |
| Application      | WB, IHC-P, ELISA                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                 |
| SWISS            | O43819                                                                            |
| Protein Weight   | 30kDa                                                                             |
| Shipping         | Ice bag                                                                           |