


Anti-Scp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22307

Item Number: CM0027031



Introduction

SCP2 Rabbit pAb validated for WB, IHC-P, ELISA in Human, Mouse. SCP2 is a peroxisomal protein that catalyzes the last step of branched-chain fatty acid beta-oxidation and mediates lipid transfer. Suitable for metabolic and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	SCP2 Rabbit pAb
Remarks/Alias	NLTP; SCOX; SCPX; SCP-2; SCP-X; NSL-TP; SCP-CHI; SCP2
Species	Human
Gene ID (Human)	6342
Gene ID	6342
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-322 of human SCP2 (NP_001007099.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P22307
Protein Weight	59kDa
Shipping	Ice bag