

Anti-Scrn3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q0Vdg4

Item Number: CM0015389

Introduction

SCRN3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat SCR3. Ideal for studying thermal nociception mechanisms.

[Learn More](#)

Parameter	Specification
Product Name	SCRN3 Rabbit pAb
Remarks/Alias	SE53; SCR3
Species	Human
Gene ID (Human)	79634
Gene ID	79634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-410 of human SCR3 (NP_078859.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q0VDG4
Protein Weight	49kDa
Shipping	Ice bag