

Anti-Sctr Rabbit Polyclonal Antibody - P47872

Item Number: CM0029023



Introduction

Rabbit polyclonal antibody against human SCTR (secretin receptor), suitable for WB and ELISA. Cross-reacts with human and mouse. Ideal for research on pH regulation, thermogenesis, and water homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SCTR Rabbit pAb
Remarks/Alias	SR; SCTR
Species	Human
GeneID (Human)	6344
GeneID	6344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 151-461 of human SCTR (P47872).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47872
Protein Weight	50kDa
Shipping	Ice bag