

Anti-Scyl2 Polyclonal Antibody For Wb And Elisa - Q6P3W7

Item Number: CM0017736

Introduction

Anti-SCYL2 rabbit polyclonal antibody validated for WB and ELISA, reacts with human and rat. SCYL2, a clathrin-associated kinase, regulates neuronal receptor trafficking and brain development.

[Learn More](#)

Parameter	Value
Product Name	SCYL2 Rabbit pAb
Remarks/Alias	AMC4; AMCNACC; CVAK104; SCYL2
Species	Human
GeneID (Human)	55681
GeneID	55681
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SCYL2 (NP_060458.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6P3W7
Protein Weight	104kDa
Shipping	Ice bag