

Anti-Sd1 Polyclonal Antibody For Wb And Elisa - B6F2D9

Item Number: CM0022120

Introduction

Rabbit polyclonal antibody against Oryza sativa SD1 (gibberellin 20-oxidase). Validated for Western blot and ELISA, ideal for studying gibberellin biosynthesis and plant growth regulation.

[Learn More](#)

Parameter	Value
Product Name	SD1 Rabbit pAb
Remarks/Alias	sd1; GA20; Sd-1; 20ox2; C20ox2; SD1-2E; GA20OX2; Os20ox2; osGA20ox2; SD1
Species	Oryza sativa
Gene ID	4325003
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of Oryza sativa SD1 (Q0JH50)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	B6F2D9
Protein Weight	43kDa
Shipping	Ice bag