

Anti-Sdc2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P34741

Item Number: CM0018860



Introduction

Anti-SDC2 rabbit pAb validated for WB, IHC-P, ELISA in human, mouse, rat. Targets syndecan-2, a heparan sulfate proteoglycan regulating dendritic arbor morphogenesis and neurogenesis.

[Learn More](#)

Parameter	Value
Product Name	SDC2 Rabbit pAb
Remarks/Alias	HSPG; CD362; HSPG1; SYND2; SDC2
Species	Human
Gene ID	6383
Immunogen	Recombinant fusion protein containing amino acids 19-144 of human SDC2 (NP_002989.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34741
Protein Weight	22kDa
Shipping	Ice bag