

Anti-Sdc4 Rabbit Polyclonal Antibody For Wb, Elisa - P31431

Item Number: CM0024147



Introduction

High-quality rabbit polyclonal antibody against human SDC4 (Syndecan-4), validated for WB and ELISA. Ideal for exosome biogenesis research and diagnostic assay development. Recognizes ~22 kDa band.

[Learn More](#)

Parameter	Value
Product Name	SDC4 Rabbit pAb
Remarks/Alias	SYND4; SDC4
Species	Human
Gene ID (Human)	6385
Gene ID	6385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-145 of human SDC4 (NP_002990.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31431
Protein Weight	22kDa
Shipping	Ice bag