

Anti-Sdha Monoclonal Antibody For Wb, If, Elisa - P31040

Item Number: CM0002697



Introduction

Rabbit recombinant monoclonal antibody targeting human SDHA, validated for Western blot, immunofluorescence, and ELISA. Detects endogenous SDHA in human, mouse, and rat samples. Suitable for mitochondrial complex II research and tumor suppression studies.

[Learn More](#)

Parameter	Value
Product Name	SDHA Rabbit mAb
Remarks/Alias	FP; PGL5; SDH1; SDH2; SDHF; CMD1GG; MC2DN1; NDAXOA; SDHA
Species	Human
Gene ID (Human)	6389
Gene ID	6389
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 550-650 of human SDHA (P31040).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31040
Protein Weight	73 kDa
Shipping	Ice bag