


## Anti-Sdha Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P31040

Item Number: CM0025912



### Introduction

Rabbit polyclonal antibody targeting human SDHA (succinate dehydrogenase complex flavoprotein subunit). Confirmed for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen covering the C-terminal region. Ideal for mitochondrial function, metabolic disease, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SDHA Rabbit pAb                                                                                                                        |
| Remarks/Alias    | FP; PGL5; SDH1; SDH2; SDHF; CMD1GG; MC2DN1; NDAXOA; SDHA                                                                               |
| Species          | Human                                                                                                                                  |
| Gene ID (Human)  | 6389                                                                                                                                   |
| Gene ID          | 6389                                                                                                                                   |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 385-664 of human SDHA (NP_004159.2). |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                      |
| SWISS            | P31040                                                                                                                                 |
| Protein Weight   | 73kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |