


Anti-Sdhaf2 Polyclonal Antibody For Wb, Elisa - Q9Nx18

Item Number: CM0022905



Introduction

SDHAF2 Rabbit pAb for WB and ELISA, targeting human SDHAF2 (mitochondrial succinate dehydrogenase assembly factor 2). Cross-reacts with mouse. Ideal for mitochondrial metabolism and TCA cycle research.

[Learn More](#)

Parameter	Value
Product Name	SDHAF2 Rabbit pAb
Remarks/Alias	PGL2; SDH5; hSDH5; C11orf79; SDHAF2
Species	Human
GeneID (Human)	54949
GeneID	54949
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-166 of human SDHAF2 (NP_060311.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NX18
Protein Weight	20kDa
Shipping	Ice bag