


Anti-Sdhb Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P21912

Item Number: CM0006542



Introduction

Recombinant rabbit monoclonal antibody targeting human SDHB. Validated for WB, IF/ICC, IHC-P, ELISA. Detects endogenous SDHB in human, mouse, rat. Ideal for mitochondrial complex II and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SDHB Rabbit mAb
Remarks/Alias	CWS2; IP; PGL4; SDH; SDH1; SDH2; SDHIP; succinate dehydrogenase complex iron sulfur subunit B
Species	Human
Gene ID (Human)	6390
Gene ID	6390
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-280 of human SDHB (P21912)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21912
Protein Weight	31kDa
Shipping	Ice bag