


Anti-Sdhb Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P21912

Item Number: CM0021426



Introduction

Rabbit polyclonal antibody targeting human SDHB, the iron-sulfur subunit of mitochondrial complex II. Validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat SDHB. Predicted molecular weight 32 kDa.

[Learn More](#)

Parameter	Value
Product Name	SDHB Rabbit pAb
Remarks/Alias	IP, SDH, CWS2, PGL4, SDH1, SDH2, SDHIP, MC2DN4, SDHB
Species	Human
GeneID (human)	6390
GeneID	6390
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-280 of human SDHB (NP_002991.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21912
Protein Weight	32kDa
Shipping	Ice bag