

Anti-Sdhc Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99643

Item Number: CM0021715



Introduction

Rabbit polyclonal antibody against human SDHC (succinate dehydrogenase complex subunit C), suitable for WB, IHC-P, IF/ICC, and ELISA.

Detects human, mouse, and rat SDHC. Target is a mitochondrial membrane protein involved in electron transport.

[Learn More](#)

Parameter	Value
Product Name	SDHC Rabbit pAb
Remarks/Alias	CYBL; PGL3; QPS1; SDH3; CYB560; SDHC
Species	Human
GeneID (Human)	6391
GeneID	6391
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SDHC (NP_002992.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99643
Protein Weight	19kDa
Shipping	Ice bag