


Anti-Sdhc Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q99643

Item Number: CM0007368



Introduction

Rabbit monoclonal antibody against human SDHC (CYBL, SDH3), validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial complex II research.

[Learn More](#)

Parameter	Value
Product Name	SDHC Rabbit mAb
Remarks/Alias	CYBL; PGL3; QPS1; SDH3; CYB560; SDHC
Species	Human
GeneID	6391
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SDHC (NP_002992.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99643
Protein Weight	19 kDa
Shipping	Ice bag