

Anti-Sdhd Rabbit Polyclonal Antibody For Wb, Elisa - O14521

Item Number: CM0023923



Introduction

Rabbit polyclonal antibody targeting SDHD (succinate dehydrogenase subunit D) for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat, enabling mitochondrial complex II research.

[Learn More](#)

Parameter	Value
Product Name	SDHD Rabbit pAb
Remarks/Alias	PGL; CBT1; CWS3; PGL1; QPs3; SDH4; cybS; CII-4; MC2DN3; SDHD
Species	Human
Gene ID (Human)	6392
Gene ID	6392
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human SDHD (NP_002993.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14521
Protein Weight	17kDa
Shipping	Ice bag