

Anti-Sdpr/Cavin-2 Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - O95810

Item Number: CM0004617

Introduction

SDPR/Cavin-2 Rabbit PolymAb® detects Human and Mouse cavin-2 (SWISS: O95810) in WB, IF-P, IHC-P, ELISA. This caveolae-associated protein is key to membrane curvature and caveolar biogenesis. Ideal for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	SDPR/Cavin-2 Rabbit PolymAb®
Remarks/Alias	SDR; SDPR; PS-p68; cavin-2
Species	Human
GeneID (Human)	8436
GeneID	8436
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 326-425 of human SDPR/Cavin-2 (NP_004648.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95810
Protein Weight	47kDa
Shipping	Ice bag