

Anti-Sdsl Polyclonal Antibody For Wb, Elisa - Q96Ga7

Item Number: CM0015935



Introduction

Anti-SDSL Rabbit Polyclonal Antibody designed for WB and ELISA applications. Targets human serine dehydratase-like protein (SWISS Q96GA7), cross-reacts with mouse, suitable for lung cancer and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	SDSL Rabbit pAb
Remarks/Alias	TDH; cSDH; SDH 2; SDS-RS1; SDSL
Species	Human
GeneID (Human)	113675
GeneID	113675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-137 of human SDSL (NP_612441.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96GA7
Protein Weight	35kDa
Shipping	Ice bag