


Anti-Sec14L2 Polyclonal Antibody For Ihc-P, Elisa - O76054

Item Number: CM0020697



Introduction

SEC14L2 (TAP/SPF) rabbit polyclonal antibody validated for IHC-P and ELISA. Cross-reacts with human, mouse, rat. For studying alpha-tocopherol/squalene transport, cholesterol biosynthesis, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	SEC14L2 Rabbit pAb
Remarks/Alias	SEC14L2; C22orf6; SPF; TAP; TAP1; SEC14-like protein 2
Species	Human
Gene ID (Human)	23541
Gene ID	23541
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-275 of human SEC14L2 (NP_203740.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76054
Protein Weight	36kDa/44kDa/46kDa
Shipping	Ice bag