


Anti-Sec14L2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O76054

Item Number: CM0016027



Introduction

Anti-SEC14L2 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat SEC14L2 (46 kDa), a lipid carrier/transcriptional activator associated with tocopherol and cholesterol.

[Learn More](#)

Parameter	Value
Product Name	SEC14L2 Rabbit pAb
Remarks/Alias	SPF; TAP; TAP1; C22orf6; SEC14L2
Species	Human
GeneID (human)	23541
GeneID	23541
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-270 of human SEC14L2 (NP_001191133.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76054
Protein Weight	46kDa
Shipping	Ice bag