


Anti-Sec23A Polyclonal Antibody For Wb, Elisa - Q15436

Item Number: CM0016749



Introduction

Rabbit polyclonal antibody against human Sec23A, validated for WB and ELISA. Cross-reacts with mouse and rat Sec23A. Suitable for COPII vesicle trafficking and ER-Golgi transport studies.

[Learn More](#)

Parameter	Value
Product Name	Sec23A Rabbit pAb
Remarks/Alias	CLSD; hSec23A; Sec23A
Species	Human
Gene ID (Human)	10484
Gene ID	10484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human Sec23A (NP_006355.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15436
Protein Weight	86kDa
Shipping	Ice bag