

Anti-Sec23A Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q15436

Item Number: CM0005761



Introduction

High-quality anti-Sec23A rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Targets human Sec23A, cross-reacts with mouse. Ideal for COPII vesicle trafficking and ER-Golgi transport studies.

[Learn More](#)

Parameter	Value
Product Name	Sec23A Rabbit mAb
Remarks/Alias	CLSD; hSec23A; Sec23A
Species	Human
GeneID (Human)	10484
GeneID	10484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human Sec23A (Q15436)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15436
Protein Weight	86kDa
Shipping	Ice bag