

Anti-Sec23B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15437

Item Number: CM0017260

Introduction

Rabbit polyclonal antibody targeting human Sec23B, a COPII coat protein essential for ER-to-Golgi trafficking. Suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat. Predicted molecular weight 86 kDa. Ideal for secretion and membrane transport studies.

[Learn More](#)

Parameter	Value
Product Name	Sec23B Rabbit pAb
Remarks/Alias	CWS7; CDAII; CDAN2; CDA-II; HEMPAS; hSec23B; Sec23B
Species	Human
GeneID (Human)	10483
GeneID	10483
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human Sec23B (NP_006354.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15437
Protein Weight	86kDa
Shipping	Ice bag