

Anti-Sec23Ip Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6Y8

Item Number: CM0015144

Introduction

Rabbit polyclonal antibody against SEC23IP (p125), validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for human, mouse, and rat samples; targets the SEC23-interacting protein involved in ER exit site organization and COPII vesicle formation.

[Learn More](#)

Parameter	Value
Product Name	SEC23IP Rabbit pAb
Remarks/Alias	P125; P125A; iPLA1A; MSTP053; iPLA1beta; SEC23IP
Species	Human
Gene ID (Human)	11196
Gene ID	11196
Immunogen	Recombinant Protein: Recombinant fusion protein containing a sequence corresponding to amino acids 600-940 of human SEC23IP (NP_009121.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6Y8
Protein Weight	111kDa
Shipping	Ice bag