


Anti-Sec24C Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P53992

Item Number: CM0015113



Introduction

Rabbit polyclonal antibody against human SEC24C, validated for Western blot, immunofluorescence, and ELISA. Targets the COPII vesicle cargo selection protein essential for ER-to-Golgi transport.

[Learn More](#)

Parameter	Value
Product Name	SEC24C Rabbit pAb
Remarks/Alias	SEC24C
Species	Human
Gene ID (Human)	9632
Gene ID	9632
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 915-1094 of human SEC24C (NP_004913.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P53992
Protein Weight	118kDa
Shipping	Ice bag