


Anti-Sec24D Monoclonal Antibody For Wb, Elisa - 094855

Item Number: CM0005948



Introduction

Rabbit monoclonal antibody targeting SEC24D, a COPII vesicle coat protein essential for ER-to-Golgi transport. Validated for WB, ELISA across human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	SEC24D Rabbit mAb
Remarks/Alias	CLCRP2; SEC24D
Species	Human
GeneID (Human)	9871
GeneID	9871
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SEC24D (NP_001304995.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94855
Protein Weight	113kDa
Shipping	Ice bag