


## Anti-Sec24D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O94855

Item Number: CM0020083



### Introduction

Rabbit polyclonal antibody against human SEC24D, validated for WB, IF/ICC, ELISA; cross-reacts with mouse, rat. Suitable for studying COPII vesicle transport and ER-Golgi cargo selection.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SEC24D Rabbit pAb                                                                                                       |
| Remarks/Alias    | CLCRP2; SEC24D                                                                                                          |
| Species          | Human                                                                                                                   |
| Gene ID          | 9871                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 830-1000 of human SEC24D (NP_001304995.1) |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, IF/ICC, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                                       |
| SWISS            | O94855                                                                                                                  |
| Protein Weight   | 113kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                 |