

Anti-Sec61A Rabbit Monoclonal Antibody For Wb, Elisa - Q9H9S3

Item Number: CM0003283

Introduction

Rabbit monoclonal antibody targeting human Sec61A (Q9H9S3) for Western blot and ELISA. Reacts with human, mouse, and rat. Essential tool for investigating protein translocation across the endoplasmic reticulum.

[Learn More](#)

Parameter	Value
Product Name	Sec61A Rabbit mAb
Remarks/Alias	Sec61 alpha-2; SEC61A2; Sec61A
Species	Human
GeneID (Human)	55176
GeneID	55176
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human Sec61A (Q9H9S3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9S3
Protein Weight	52kDa
Shipping	Ice bag