


Anti-Sec61A1 Polyclonal Antibody For Wb, Elisa - P61619

Item Number: CM0017967



Introduction

SEC61A1 Rabbit pAb validated for WB and ELISA. Recognizes human SEC61A1, key for ER protein translocation and nephrogenesis, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	SEC61A1 Rabbit pAb
Remarks/Alias	HNFJ4; SEC61; ADTKD5; HSEC61; SEC61A; SEC61A1
Species	Human
GeneID (Human)	29927
GeneID	29927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-420 of human SEC61A1 (NP_037468.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P61619
Protein Weight	52kDa
Shipping	Ice bag