


Anti-Sec61A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P61619

Item Number: CM0015765



Introduction

Rabbit polyclonal antibody against human SEC61A1 (P61619), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Used for studying ER protein translocation and calcium homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SEC61A1 Rabbit pAb
Remarks/Alias	HNFJ4; SEC61; ADTKD5; HSEC61; SEC61A; SEC61A1
Species	Human
Gene ID	29927
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-420 of human SEC61A1 (NP_037468.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61619
Protein Weight	52kDa
Shipping	Ice bag