


Anti-Sec61B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P60468

Item Number: CM0022067



Introduction

High-quality rabbit polyclonal antibody targeting human SEC61B, suitable for Western Blot, IF/ICC, and ELISA. Detects human, mouse, and rat SEC61B. Ideal for studying ER protein translocation and membrane insertion, ensuring reliable performance for your protein transport research.

[Learn More](#)

Parameter	Value
Product Name	SEC61B Rabbit pAb
Remarks/Alias	SEC61B; Sec61 translocon beta subunit
Species	Human
Gene ID (Human)	10952
Gene ID	10952
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-96 of human SEC61B (NP_006799.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60468
Protein Weight	10kDa
Shipping	Ice bag