


Anti-Sec62 Polyclonal Antibody For Wb, Elisa - Q99442

Item Number: CM0019660



Introduction

High-quality SEC62 Rabbit Polyclonal Antibody validated for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat samples. Ideal for research on ER translocation and protein transport mechanisms.

[Learn More](#)

Parameter	Specification
Product Name	SEC62 Rabbit pAb
Remarks/Alias	HTP1; TP-1; Dtrp1; TLOC1; SEC62
Species	Human
GeneID	7095
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 258-399 of human SEC62 (NP_003253.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99442
Protein Weight	46kDa
Shipping	Ice bag