


Anti-Sec62 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q99442

Item Number: CM0006896



Introduction

SEC62 rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Targets human, mouse, rat translocation protein SEC62 (~46kDa). Ideal for ER protein transport research. UniProt Q99442.

[Learn More](#)

Parameter	Value
Product Name	SEC62 Rabbit mAb
Remarks/Alias	HTP1; TP-1; Dtrp1; TLOC1; SEC62
Species	Human
Gene ID (Human)	7095
Gene ID	7095
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 214-314 of human SEC62 (NP_003253.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99442
Protein Weight	46kDa
Shipping	Ice bag