

Anti-Sec63 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ugp8

Item Number: CM0019405

Introduction

Rabbit polyclonal antibody against human SEC63 (ERdj2, DNAJC23), validated for WB and ELISA. Cross-reacts with human, mouse, rat. Targets the ER membrane translocon component for protein trafficking research.

[Learn More](#)

Parameter	Value
Product Name	SEC63 Rabbit pAb
Remarks/Alias	ERdj2; PCLD2; SEC63L; DNAJC23; PRO2507; SEC63
Species	Human
GeneID (Human)	11231
GeneID	11231
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 595-720 of human SEC63 (NP_009145.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UGP8
Protein Weight	88kDa
Shipping	Ice bag